



SNRNP40 Antibody, HRP conjugated

Product Code	CSB-PA822196LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DI7
Immunogen	Recombinant Human U5 small nuclear ribonucleoprotein 40 kDa protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U5 small nuclear ribonucleoprotein 40 kDa protein (U5 snRNP 40 kDa protein) (U5-40K) (38 kDa-splicing factor) (Prp8-binding protein) (hPRP8BP) (U5 snRNP-specific 40 kDa protein) (WD repeat-containing protein 57), SNRNP40, PRP8BP SFP38 WDR57
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRNP40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.