




FAHD2A Antibody

Product Code	CSB-PA822218LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GK7
Immunogen	Recombinant Human Fumarylacetoacetate hydrolase domain-containing protein 2A protein (4-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fumarylacetoacetate hydrolase domain-containing protein 2A (EC 3.-.-.-), FAHD2A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAHD2A

Image



Western blot

All lanes: FAHD2A antibody at 0.5µg/ml

Lane 1: K562 whole cell lysate

Lane 2: U87 whole cell lysate


Lane 3: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA822218LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.