



FAHD2A Antibody, Biotin conjugated

Product Code	CSB-PA822218LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GK7
Immunogen	Recombinant Human Fumarylacetoacetate hydrolase domain-containing protein 2A protein (4-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fumarylacetoacetate hydrolase domain-containing protein 2A (EC 3.-.-.-), FAHD2A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.