



# ZDHHC12 Antibody, Biotin conjugated

|                     |                                                                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA822219LD01HU                                                                                                                                        |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                             |
| Uniprot No.         | Q96GR4                                                                                                                                                    |
| Immunogen           | Recombinant Human Probable palmitoyltransferase ZDHHC12 protein (121-197AA)                                                                               |
| Raised In           | Rabbit                                                                                                                                                    |
| Species Reactivity  | Human                                                                                                                                                     |
| Tested Applications | ELISA                                                                                                                                                     |
| Form                | Liquid                                                                                                                                                    |
| Conjugate           | Biotin                                                                                                                                                    |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                          |
| Purification Method | Antigen Affinity Purified                                                                                                                                 |
| Isotype             | IgG                                                                                                                                                       |
| Clonality           | Polyclonal                                                                                                                                                |
| Alias               | Probable palmitoyltransferase ZDHHC12 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 12) (DHHC-12) (Zinc finger protein 400), ZDHHC12, ZNF400 |
| Immunogen Species   | Homo sapiens (Human)                                                                                                                                      |
| Research Area       | Neuroscience                                                                                                                                              |