



HHIPL1 Antibody

Product Code	CSB-PA822236LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JK4
Immunogen	Recombinant Human HHIP-like protein 1 protein (401-608AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HHIP-like protein 1, HHIPL1, HHIP2 KIAA1822
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HHIPL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.