



MS4A4A Antibody, HRP conjugated

Product Code	CSB-PA822238LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JQ5
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 4A protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 4A (CD20 antigen-like 1) (Four-span transmembrane protein 1), MS4A4A, 4SPAN1 CD20L1 MS4A4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.