



MS4A4A Antibody, HRP conjugated

Product Code	CSB-PA822238LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JQ5
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 4A protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 4A (CD20 antigen-like 1) (Four-span transmembrane protein 1), MS4A4A, 4SPAN1 CD20L1 MS4A4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A4A