



OPALIN Antibody, HRP conjugated

Product Code	CSB-PA822279LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PE5
Immunogen	Recombinant Human Opalin protein (51-141AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opalin (Oligodendrocytic myelin paranodal and inner loop protein) (Transmembrane protein 10), OPALIN, HTMP10 TMEM10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OPALIN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.