



TTC5 Antibody, HRP conjugated

Product Code	CSB-PA822676LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0Z6
Immunogen	Recombinant Human Tetratricopeptide repeat protein 5 protein (201-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 5 (TPR repeat protein 5) (Stress-responsive activator of p300) (Strap), TTC5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.