




C4orf33 Antibody

Product Code	CSB-PA822682LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1A6
Immunogen	Recombinant Human UPF0462 protein C4orf33 protein (1-199AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0462 protein C4orf33, C4orf33
Immunogen Species	Homo sapiens (Human)
Research Area	Others

Image



Western blot
 All lanes: C4orf33 antibody at 4µg/ml
 Lane 1: Mouse kidney tissue
 Lane 2: Mouse liver tissue
 Lane 3: Mouse lung tissue
 Lane 4: Mouse heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 24 kDa
 Observed band size: 24 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.