





BATF2 Antibody

Product Code	CSB-PA822685LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1L9
Immunogen	Recombinant Human Basic leucine zipper transcriptional factor ATF-like 2 protein (1-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Basic leucine zipper transcriptional factor ATF-like 2 (B-ATF-2) (Suppressor of AP-1 regulated by IFN) (SARI), BATF2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	BATF2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA822685LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA822685LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.