



# OR2L13 Antibody, HRP conjugated

Product Code	CSB-PA822694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N349
Immunogen	Recombinant Human Olfactory receptor 2L13 protein (158-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2L13 (Olfactory receptor 2L14), OR2L13, OR2L14
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	OR2L13