



OR2L13 Antibody, Biotin conjugated

Product Code	CSB-PA822694LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N349
Immunogen	Recombinant Human Olfactory receptor 2L13 protein (158-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2L13 (Olfactory receptor 2L14), OR2L13, OR2L14
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	OR2L13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.