



PPM1K Antibody, HRP conjugated

Product Code	CSB-PA822697LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3J5
Immunogen	Recombinant Human Protein phosphatase 1K, mitochondrial protein (30-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1K, mitochondrial (EC 3.1.3.16) (PP2C domain-containing protein phosphatase 1K) (PP2C-like mitochondrial protein) (PP2C-type mitochondrial phosphoprotein phosphatase) (PTMP) (Protein phosphatase 2C isoform kappa) (PP2C-kappa), PPM1K, PP2CM
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PPM1K
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.