



# PPM1K Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA822697LC01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8N3J5  |
| <b>Immunogen</b>           | Recombinant Human Protein phosphatase 1K, mitochondrial protein (30-372AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Protein phosphatase 1K, mitochondrial (EC 3.1.3.16) (PP2C domain-containing protein phosphatase 1K) (PP2C-like mitochondrial protein) (PP2C-type mitochondrial phosphoprotein phosphatase) (PTMP) (Protein phosphatase 2C isoform kappa) (PP2C-kappa), PPM1K, PP2CM |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Tags & Cell Markers   |
| <b>Target Names</b>        | PPM1K   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |