



SPEM1 Antibody, Biotin conjugated

Product Code	CSB-PA822704LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4L4
Immunogen	Recombinant Human Spermatid maturation protein 1 protein (50-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatid maturation protein 1, SPEM1, C17orf83
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	SPEM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.