



FAM160A2 Antibody, Biotin conjugated

Product Code	CSB-PA822711LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N612
Immunogen	Recombinant Human FTS and Hook-interacting protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FTS and Hook-interacting protein (FHIP) (Protein FAM160A2), FAM160A2, C11orf56 KIAA1759
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM160A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.