



# FAM160A2 Antibody, Biotin conjugated

Product Code	CSB-PA822711LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N612
Immunogen	Recombinant Human FTS and Hook-interacting protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FTS and Hook-interacting protein (FHIP) (Protein FAM160A2), FAM160A2, C11orf56 KIAA1759
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM160A2