



# PCDH20 Antibody, HRP conjugated

Product Code	CSB-PA822719LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6Y1
Immunogen	Recombinant Human Protocadherin-20 protein (61-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin-20 (Protocadherin-13), PCDH20, PCDH13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PCDH20