



TERB1 Antibody, HRP conjugated

Product Code	CSB-PA822742LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NA31
Immunogen	Recombinant Human Telomere repeats-binding bouquet formation protein 1 protein (301-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Telomere repeats-binding bouquet formation protein 1 (Coiled-coil domain-containing protein 79), TERB1, CCDC79
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TERB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.