



C12orf50 Antibody

Product Code	CSB-PA822743LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NA57
Immunogen	Recombinant Human Uncharacterized protein C12orf50 protein (201-414AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C12orf50, C12orf50
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C12orf50
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.