



NSUN7 Antibody, HRP conjugated

Product Code	CSB-PA822769LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE18
Immunogen	Recombinant Human Putative methyltransferase NSUN7 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative methyltransferase NSUN7 (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 7), NSUN7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NSUN7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.