



TCP10L Antibody, Biotin conjugated

Product Code	CSB-PA823465LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDR4
Immunogen	Recombinant Human T-complex protein 10A homolog 2 protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-complex protein 10A homolog 2 (T-complex protein 10A-2) (TCP10A-2) (TCP10-like), TCP10L
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TCP10L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.