



SLC22A17 Antibody, HRP conjugated

Product Code	CSB-PA823881LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUG5
Immunogen	Recombinant Human Solute carrier family 22 member 17 protein (1-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 17 (24p3 receptor) (24p3R) (Brain-type organic cation transporter) (Lipocalin-2 receptor) (Neutrophil gelatinase-associated lipocalin receptor) (NgalR), SLC22A17, BOCT BOIT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A17