



SLC22A17 Antibody, Biotin conjugated

Product Code	CSB-PA823881LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUG5
Immunogen	Recombinant Human Solute carrier family 22 member 17 protein (1-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 17 (24p3 receptor) (24p3R) (Brain-type organic cation transporter) (Lipocalin-2 receptor) (Neutrophil gelatinase-associated lipocalin receptor) (NgalR), SLC22A17, BOCT BOIT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.