



# PTCD2 Antibody, Biotin conjugated

Product Code	CSB-PA823887LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV60
Immunogen	Recombinant Human Pentatricopeptide repeat-containing protein 2, mitochondrial protein (153-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pentatricopeptide repeat-containing protein 2, mitochondrial, PTCD2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTCD2