




C4orf3 Antibody

Product Code	CSB-PA823892LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVX3
Immunogen	Recombinant Human Uncharacterized protein C4orf3 protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C4orf3 (Hepatitis C virus F protein-transactivated protein 1) (HCV F-transactivated protein 1), C4orf3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology

Image



Western Blot

All lanes: C4orf3 antibody at 6µg/ml

Lane 1: 293T whole cell lysate

Lane 2: Hela whole cell lysate

Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 8, 23 kDa

Observed band size: 8 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.