



Serpina3n Antibody, HRP conjugated

Product Code	CSB-PA835578LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91WP6
Immunogen	Recombinant Mouse Serine protease inhibitor A3N protein (21-418AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor A3N (Serpina A3N), Serpina3n, Spi2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Serpina3n