



Serpina3n Antibody, Biotin conjugated

Product Code	CSB-PA835578LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91WP6
Immunogen	Recombinant Mouse Serine protease inhibitor A3N protein (21-418AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor A3N (Serpin A3N), Serpina3n, Spi2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Serpina3n
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.