



SUSD6 Antibody, FITC conjugated

Product Code	CSB-PA835675LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92537
Immunogen	Recombinant Human Sushi domain-containing protein 6 protein (144-303AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sushi domain-containing protein 6 (Drug-activated gene overexpressed protein), SUSD6, DRAGO KIAA0247
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SUSD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.