



RFK Antibody, Biotin conjugated

Product Code	CSB-PA836167LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969G6
Immunogen	Recombinant Human Riboflavin kinase protein (1-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Riboflavin kinase (EC 2.7.1.26) (ATP:riboflavin 5'-phosphotransferase) (Flavokinase), RFK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RFK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.