



PNPLA2 Antibody, Biotin conjugated

Product Code	CSB-PA836180LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AD5
Immunogen	Recombinant Human Patatin-like phospholipase domain-containing protein 2 protein (1-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Patatin-like phospholipase domain-containing protein 2 (EC 3.1.1.3) (Adipose triglyceride lipase) (Calcium-independent phospholipase A2) (Desnutrin) (IPLA2-zeta) (Pigment epithelium-derived factor) (TTS2.2) (Transport-secretion protein 2) (TTS2), PNPLA2, ATGL
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PNPLA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.