



RGMA Antibody, Biotin conjugated

Product Code	CSB-PA836186LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96B86
Immunogen	Recombinant Human Repulsive guidance molecule A protein (153-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Repulsive guidance molecule A (RGM domain family member A), RGMA, RGM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RGMA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.