




TRIM43 Antibody

Product Code	CSB-PA836189LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BQ3
Immunogen	Recombinant Human Tripartite motif-containing protein 43 protein (43-240AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 43, TRIM43
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM43

Image



Western Blot

Positive WB detected in: Mouse brain tissue

All lanes: TRIM43 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa