



FAM167A Antibody, HRP conjugated

Product Code	CSB-PA836244LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KS9
Immunogen	Recombinant Human Protein FAM167A protein (1-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM167A, FAM167A, C8orf13
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM167A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.