



C8orf59 Antibody

Product Code	CSB-PA836635LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0T1
Immunogen	Recombinant Human Uncharacterized protein C8orf59 protein (1-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C8orf59, C8orf59
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C8orf59