



C8orf59 Antibody

Product Code	CSB-PA836635LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0T1
Immunogen	Recombinant Human Uncharacterized protein C8orf59 protein (1-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein C8orf59, C8orf59
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C8orf59
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.