




TTC9C Antibody

Product Code	CSB-PA836671LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5M4
Immunogen	Recombinant Human Tetratricopeptide repeat protein 9C protein (1-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 9C (TPR repeat protein 9C), TTC9C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTC9C

Image



Western blot


All lanes: TTC9C antibody at 12µg/ml + THP-1 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 21, 25 kDa

Observed band size: 21 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA836671LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.