



TTC9C Antibody, HRP conjugated

Product Code	CSB-PA836671LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5M4
Immunogen	Recombinant Human Tetratricopeptide repeat protein 9C protein (1-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 9C (TPR repeat protein 9C), TTC9C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTC9C