



L3MBTL4 Antibody, Biotin conjugated

Product Code	CSB-PA836703LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NA19
Immunogen	Recombinant Human Lethal(3)malignant brain tumor-like protein 4 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lethal(3)malignant brain tumor-like protein 4 (H-I(3)mbt-like protein 4) (L(3)mbt-like protein 4) (L3mbt-like 4), L3MBTL4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	L3MBTL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.