



ST6GALNAC3 Antibody, Biotin conjugated

Product Code	CSB-PA836724LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDV1
Immunogen	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 protein (40-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 (EC 2.4.99.-) (GalNAc alpha-2,6-sialyltransferase III) (ST6GalNAc III) (ST6GalNAcIII) (STY) (Sialyltransferase 7C) (SIAT7-C), ST6GALNAC3, SIAT7C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ST6GALNAC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.