




PPTC7 Antibody

Product Code	CSB-PA836783LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI37
Immunogen	Recombinant Human Protein phosphatase PTC7 homolog protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase PTC7 homolog (EC 3.1.3.16) (T-cell activation protein phosphatase 2C) (TA-PP2C) (T-cell activation protein phosphatase 2C-like), PPTC7, TAPP2C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPTC7

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA836783LA01HU at dilution of 1:100