



CXorf40A Antibody, FITC conjugated

Product Code	CSB-PA837446LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE69
Immunogen	Recombinant Human Protein CXorf40A protein (21-158AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein CXorf40A (Endothelial-overexpressed lipopolysaccharide-associated factor 1), CXorf40A, CXorf40 EOLA1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CXorf40A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.