



RNF128 Antibody, Biotin conjugated

Product Code	CSB-PA837448LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEB7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF128 protein (229-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF128 (EC 2.3.2.27) (Gene related to anergy in lymphocytes protein) (GRAIL) (RING finger protein 128) (RING-type E3 ubiquitin transferase RNF128), RNF128
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF128
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.