



RNF128 Antibody, Biotin conjugated

Product Code	CSB-PA837448LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEB7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF128 protein (229-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF128 (EC 2.3.2.27) (Gene related to anergy in lymphocytes protein) (GRAIL) (RING finger protein 128) (RING-type E3 ubiquitin transferase RNF128), RNF128
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF128