



FATE1 Antibody, Biotin conjugated

Product Code	CSB-PA839274LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969F0
Immunogen	Recombinant Human Fetal and adult testis-expressed transcript protein (1-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fetal and adult testis-expressed transcript protein (Cancer/testis antigen 43) (CT43) (Tumor antigen BJ-HCC-2), FATE1, FATE
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.