



XAGE2 Antibody, Biotin conjugated

Product Code	CSB-PA839331LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GT9
Immunogen	Recombinant Human X antigen family member 2 protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	X antigen family member 2 (XAGE-2) (Cancer/testis antigen 12.2) (CT12.2) (G antigen family D member 3), XAGE2, GAGED3 XAGE2B
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	XAGE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.