




COG3 Antibody

Product Code	CSB-PA839348LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JB2
Immunogen	Recombinant Human Conserved oligomeric Golgi complex subunit 3 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Conserved oligomeric Golgi complex subunit 3 (COG complex subunit 3) (Component of oligomeric Golgi complex 3) (Vesicle-docking protein SEC34 homolog) (p94), COG3, SEC34
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	COG3

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA839348LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.