



LLCFC1 Antibody, HRP conjugated

Product Code	CSB-PA839358LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96L11
Immunogen	Recombinant Human LLLL and CFNLAS motif-containing protein 1 protein (29-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LLLL and CFNLAS motif-containing protein 1 (MSSP-binding protein CTM-1), LLCFC1, C7orf34
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.