



MRGPRX2 Antibody, HRP conjugated

Product Code	CSB-PA839361LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB1
Immunogen	Recombinant Human Mas-related G-protein coupled receptor member X2 protein (286-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mas-related G-protein coupled receptor member X2, MRGPRX2, MRGX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MRGPRX2